

## SUPPLEMENTARY TABLE

**Supplementary Table 1. Primers for real-time PCR in this study.**

Gene	Forward primers (5'→3')	Reverse primers (5'→3')
RORC	GAGAAGGACAGGGAGCCAAG	GCGGAAGAAGCCCTTGAC
lncRNA SPIRE1	TTTCCTCAACTGGAGCTGGT	TCGTGCTTCCTGTATTGCAGA
FOXP3	CCCACACTGCCCTAGTCAT	TCCACCGTTGAGAGCTGGT
miR-181a-5p	AACAUUCAACGCUGUCGGUGAGU	Universal PCR Reverse Primer (#B532451; Sangon Biotech Co., Ltd.)
PRLR	CACTCCTTCTCCCTCTTTCTGG	ATAGGAGAGTTCTTTAGTTTTGCCA
IGF2BP2	AGTGGAATTGCATGGGAAAATCA	CAACGGCGGTTTCTGTGTC
KLF6	GGCAACAGACCTGCCTAGAG	CTCCCGAGCCAGAATGATTTT
SIRT1	TAGCCTTGTCAGATAAGGAAGGA	ACAGCTTCACAGTCAACTTTGT
EN2	CCGGCGTGGGTCTACTGTA	CCTCTTTGTTCCGGGTTCTTCTT
GAPDH	AACTTTGGCATTGTGGAAGG	ATTGGGGGTAGGAACA
U6	CTCGCTTCGGCAGCACA	ACGCTTCACGAATTTGCGT